


Overview of Instrument



You will find and leave the instrument (if not otherwise agreed)

- In daily shut down mode
- With a sheath container filled for at least 2.5h sorting
- An extra sheath container, filled, but not coupled to the instrument
- Nozzle, side stream illuminator, deflection plates clean
- Waterbath off
- Nozzle, clean, in its plastic box, in cabinet
- Table and instrument clean
- check the log book for sheath info

Overview of Instrument



The 100 µm Nozzle

The Cleaned & dry Nozzle placed in its plastic container. Container placed in top compartment of cabinet. 2 spare nozzle available in compartment's back. BioVis' standard nozzle found in front of compartment.

Nozzle with rubber O-ring

Nozzle holder, contains electronics



The Sample chamber

Do not open the transparent door of the sample chamber when the sample chamber (marked by 2 red dots) is moving up or down



Opening the door, whilst the moving, will generate a "full stop error" and the whole repertoire of debubbling, QC and drop delay is needed thereafter (15min)

You can open the chamber door when the green light is NOT blinking.


To put more attention the the green light bar, three dots were added, 2 red flanking one green. The green light will also blink when sample is been taken up. The door can be opened as well under that condition.

Sample chamber can be up, for exchanging tubes, or down, sealing the chamber for Run/Record/Sort.



Make sure the tube is upright not tilted.

The tube holder is hold in place by a magnet (and can be taken out). Do not load the sample with a tilted tube
Make also sure your tube does not have lid



Overview of Software

The lower bar of the Software gives important information

Instrument connected and ready

Info about a log file

opens stream pop-up window

Connected Ready [2020-06-19 10:52:41] Acquisition Completed Sorter Status

Stream mode (default/ straight down)

LED Light switch for Sample chamber and Sort vessel chamber

Indicator of Shut down fluid filling

Indicator of waste filling

Indicator of available sheath fluid

Menus of the software

File	Cytometer	QC/Standardization	Sorting	Settings	Advanced
New Experiment...	System Startup...	Start QC/Standardization	Sort Calibration	Set Channel	Delay Setting
New Experiment from Template...	System Shutdown...	Close QC/Standardization	Manual Drop Delay	Set Label	Laser Setting
New Experiment from FCS...	Long Term Shutdown...		Manual Side Stream Calibration	Set Customized Parameters	Event Rate Setting
New Compensation	Turn On		Sorter Status	Compensation Matrix	Sort Guardband Setting
Open Experiment	Turn Off		Side Stream Monitor	Compensation Library	Default Amplitude Setting
Open Compensation	Backflush		Stream Mode Switch	Events Display Setting	Collection Device Library
Convert CytExpert Experiment	Sheath Filter De-bubble		Reset Cyclone	Language Setting	Sort Mode Library
Save	Flow Cell De-bubble		Tube Position Setting	Options	Maintenance
Save As	Daily Clean...				Bubble Detector Calibration
Save As Template	Flow Cell Clean...				Sheath Tank Scale Reset
Import FCS File	Aseptic Clean...				
Export FCS File	Reset Sample Chamber				
Recent	Acq.Setting...				
Recent Template	Acq.Setting Catalog...				
Recent Compensation	Detector Configuration				
Exit	Sample Flow Parameter Adjustment				
	Cytometer Information...				

The CytExpert Software will not be explained here anyfurther
Please check the manual for CytExpert to get to know how to Set up an experiment

Start-up

10 min


Prerequisite:

- A) Instrument has two Sheath tanks, one coupled to Instrument , one refilled for later use
- B) Waste is empty
- C) Shut down fluid is not empty
- D) Deflection plates ad illuminator are clean

Decision point: Decide whether you want a plate or tube sort .

- This can be done later, but will need a Standby and debubble.
- The lower bar of the software gives you important information, in wich sort mode the doftware is.
- Please cross check with the actual position of the waste drawer.
- Changes of sort mode must be done in software AND the waste drawer
- Sortmode change Software: Sorting>Stream mode switch (pop-up window instruction)
- Sortmode change Instrument: STRAIGHT DOWN (LEFT), DEFAULT (RIGHT)

Straight Down Sort	Default Sort
<ul style="list-style-type: none"> • Straight down modus • Centre stream is used for sort • Sorted drops lets are NOT charged • 96 deep well and 384 only supported by straight down modus 	<ul style="list-style-type: none"> • Side stream modus • Side stream droplets will be charged • Up to four side streams are used for sort • Tube sorts • plate sort up to 96 well plate



Following the pop up window instructions carefully drag out the waste drawer (magnet holding) And place it to the left or right for straight down or default sort set up

Having a mismatch between Software and Hardware
Will cause flooding (the instrument) and erroneus sort


Start-up

10 min

Prerequisite:

- A) Instrument has two Sheath tanks, one coupled to Instrument , one refilled for later use
- B) Waste is empty
- C) Shut down fluid is not empty
- D) Deflection plates and illuminator are clean

- 1) Switch ON: Computer and Waterbath
- 2) CytExpert > Cytometer > Turn ON
- 3) Have the sorter status open to see stream
- 4) >Cytometer>Start up (Follow instructions of pop up window)
- 5) After start up : do twice a >Cytometer> Sheath filter de-bubble
- 6) Followed by check of any air bubbles in filter
- 7) Followed by twice cytomter >Flow Cell Debubble
- 8) Check stream, satellite droplets* fuse to droplets (often top)
- 9) QC> QC/Standardization >Start QC > Run QC /drop delay *
- 10) Drop delay separately done via >Sorting>Sort Calibration
- 11) You still can do the debubbling from here
- 12) Run Droplet break off point
- 13) Mark the green line of break off point with red line to indicate position changes/retain




QC / drop delay is not explained in this manual, but will be shown in hands-on session. → Automatic /Manual way

Have the Sorter Status window always open

- **Check that Auto Maintain is ON**
- Monitor that break off point is overlaying with the green line

In case the stream must be stopped and restarted

- Check that break off point is back at position marked by green line
- IF NOT: RUN **>Sorting>Sort Calibration (Pop up window instructions)**
- Auto-Maintain



Sorter Status

General

QC Last Run: 2020-12-31 09:55:18 Result: QC Passed.

Sort Calibration Last Run: 2020-12-31 10:01:25 Result: Passed

Drop Delay: 33.50

Instrument

Power USB: Power On

Camera Status: Connect

Sort Protection Door: Closed

Sample Station Door: Closed

Side Stream Illuminator: On

Drop Delay Mode: Auto

Side Stream Mode: Default

Sheath Pressure: 15.00 PSI

Sample Pressure: 15.57 PSI

Sheath Temperature: 23.1°C/73.5°F

Frequency: 34,650 Hz

Amplitude: 92.93 V

Auto Maintain Status: Auto Maintaining

⊞ Fewer details

Setting up Experiment and Sorting

Explanation of Sort modes

Single Mode - when having only one event per drop is the most important aspect of the sort.

Purity Mode - when the purity of the sort is the most important.

Purity1-2 Mode - when the purity of the sort is the most important. Purity 1-2 Mode is more inclusive than Purity Mode.

Enrich Mode - when recovery is the most important aspect of the sort. With Enrich, all positive events (an event that falls within the sort logic) are sorted

NOTE

Beckman Coulter recommends the Purity 1-2 Mode for sorting macro-particles.

Image to the right illustrates actually the adjacent stream BEFORE the sort droplets appear. Direction of Flow is indicated. **Target** population is red dot, **contamination** is blue. If an event is not in the centre of the droplet, it might end up in the preceding or trailing droplet. Depending on sort mode, such "just in case" droplet (preceding or trailing) is sorted as well.

Single mode: no contaminant gets sorted. Strictly one drop let is sorted, target needs to be droplet-centered, and only ONE target event per droplet. Is accepted. contamination should not end up in such "target"-droplet



Purity mode: Kind of a "wider Single mode". Target can be 1 and more per droplet, no contaminant will be sorted.

Purity 1-2 mode: No contaminant gets sorted. Preceding/trailing droplets will get sorted.

Enrich mode: Contaminant will end up in sort. When ever a target can get sorted it will get sorted, irresoective whether or not a contaminant is present.

Each sorted droplet will be charged in DEFAULT mode. In STRAIGHT down mode all non-sorted droplets are charged.

BioVis is using guard band 25%, for best purity. 15% is also available (0.1% error)



Setting up Experiment and Sorting

CytExpert is considered self-explanatory software, so only brief info will be given here.

The CytExpert Software will not be explained here any further. Please check the manual for CytExpert to get to know how to set up an experiment.

1. Set up a Compensation experiment (or have a matrix ready to apply to actual Experiment)
2. Set up the experiment with plots and gates
3. Open the "Acquisition set up" window and load the default 100 file, to have all values on 100
4. Tubes (in software) with data cannot be used setting up sort
5. Have an empty tube to apply sort set up

Sorting gates can only be **8 steps deep** into the gating hierarchy. **Autogates cannot be used** for Sort gates.





Tube or plate sort

Tube sort layout

Check that droplets enter the sort tubes in case alignment is needed:

Sorting > Tube Position Setting


More alignment options are manually over the outer tube inserts.



Sort Stream Preference

Place most precious and/or rare events in outer stream (L2 or R2) for maximum purity.

Place sorting events of high abundance or larger size (above 15µm) to inner streams (L1 or R1)



Setting up Experiment and Sorting



Tube or plate sort

Plate sort layout

Plate alignment is needed to ensure that droplets hit the wells
>Advanced > Collection Device Library.
 Then select **Calibrate** from the Collection Device Library window. Follow instructions of Pop-up window.








Plate sort layout is self explanatory. Hover over icons And drop down menu to Get more information One can group sort set ups, to Have them reused.



In case one plate holds sorts of rare events and abundant – make sure to sort the abundant before the rare events, as more sample will be taken up to "find" a rare event

Index Sorting selected: entire index sort data (max 30min) is recorded. Lets you also monitor which cells have ended up in the well (contaminant included in case...) Special gating are applied over the wells for that purpose.
 Use the drawing **Index sorting tool**



Exporting Data and Sort reports

Export your data different options

- >File>Export FCS file
- Data to pdf
- Sort reports
- Statistic reports



Keeping it clean

For proper sorting the instrument must be able to take up sample, have a stable stream (with no satellite droplets) and must be able to properly charge droplets .

- Clogs and bubbles → **Debubble function** , Cytometer> Sheath Filter / Flow Cell debubble
- Dirty deflection plates → open sort chamber remove illuminator and plates
- Dirty camera windows → open sort chamber remove illuminator
- Clogged sample line I → low or no events, BACKFLUSH,
- Clogged Sample line II → Run Cleanz full speed,

When sorting is done

1. RUN CLEANZ in your experiment, followed by MiliQ water : No events of your sample must appear
2. Never leave the instrument without DAILY CLEAN, Cytometer > Daily Clean **3 min each**


Shut Down

1. RUN **DAILY CLEAN** >Cytometer>Daily Clean run , minimum **3 min** for each step
2. RUN **SHUT DOWN**, follow shut down procedure, Clean all necessary part
3. When done, **Switch OFF instrument** , Cytometer> Switch Off and **Waterbath**

Nozzle: use the clams to release

Stream Illuminator: Drag out

Deflection plates: Drag out



Leaving the instrument to next customer on same day

Instrument is clean and on Standby

Some extra information

When handling

The instrument

- Don't open the sample door chamber when Sample chamber unit is moving

This cause way to much trouble ;)*

the sheath container

- make sure the lid is tight
- filling it up only until upper welding line

Pop-up windows

- Certain pop-up block the software
- Press hide button on the pop up to access software
- Find pop up window info in lower left of software

* Sample Chamber door

Do not open it if sample chamber, which seals the sample holder is in move
Opening the door, whilst the moving, puts instrument into Standby
Usually the full repertoire of Debubbling, QC and drop delay is needed thereafter

Once system is ready, always make sure the Window "Sorter Status" shows That "Auto Maintain" is on to monitor the stream automatic

BACKFLUSH : Cytometer> Backflush


Clean the nozzle by

- Spray with EtOH
- Dab dry
- Use a syringe with MiliQ
- Place syringe tip on top of Nozzle O ring
- Press water through nozzle
- Stream visible
- Dab dry

After start ups or Stream issues

- Check for bubbles in sheath filter
- Cytometer > Sheath filter debubbling 2x
- Cytometer > Flow Cell debubbling 2x
- Proceed with drop delay >Sorting>Sort Calibration

IMPORTANT to have activated



Sorter Status	
General	
QC Last Run:	2020-12-31 09:55:18 Result: QC Passed.
Sort Calibration Last Run:	2020-12-31 10:01:25 Result: Passed
Drop Delay:	33.50
Instrument	
Power USB:	Power On
Camera Status:	Connect
Sort Protection Door:	Closed
Sample Station Door:	Closed
Side Stream Illuminator:	On
Drop Delay Mode:	Auto
Side Stream Mode:	Default
Sheath Pressure:	15.00 PSI
Sample Pressure:	15.57 PSI
Sheath Temperature:	23.1°C/73.5°F
Frequency:	34,650 Hz
Amplitude:	92.93 V
Auto Maintain Status:	● Auto Maintaining Turn Off

⊖ Fewer details