





LEICA Stellaris 5 Confocal microscope Manual/Quick guide




Starting the microscope

1. **Turn ON** the computer, wait till Windows is booted
2. **Log in** with your UU AKKA account
3. **Turn ON** the scanner box switches below the table:
 - i. Power button ON
 - ii. Laser button ON
 - iii. Turn the laser interlock key to vertical position (emission light ON)




4. **Wait** approx. 3 minutes
5. **Start** the LAS X software and start "machine" configuration and "DMI8" microscope:




The microscope stand

Visualizing the sample through the eyepiece using fluorescent light (IL)
If you don't see any signal through the eyepiece, make sure that you are using the following settings on the TFT screen:




Visualizing the sample through the eyepiece using brightfield light (TL)
If you don't see any signal through the eyepiece, make sure that you are using the following settings on the TFT screen:





Follow these steps to have light through the eyepiece:

- Select the desired objective and use the correct immersion medium for it. On the TFT screen the objective buttons blink when there's a change in the immersion medium. Press them again after changed the immersion medium.
- Port: Eyepieces 100%
- Magn.-Changer: 1




To change the intensity for the TL or IL:



Using the LAS X software

In LAS X, hover the mouse on any button to get info about it.

Configuration tab



Under the configuration tab, several options can be changed and lasers can be turned ON and OFF.

- In the Laser config menu: laser can be turned ON or OFF. We have two lasers, a fixed 405nm and a tunable White light laser (WLL) 485-790nm.
- In the Hardware menu: dynamic range can be set to 16-bit (the default is 8-bit)
- In the USB Panel, the sensitivity of the knobs can be changed

Acquire tab, confocal imaging

This is the most used tab, here the imaging setup can be configured and imaging can be started.

▶; ▾ : Expand or collapse a menu.

Info and help menu.
Use it when you need!

Acquisition mode:
Type of imaging,
2D, 3D, time-lapse,
spectral scanning

Navigator. Quick spiral
scanning, tile scanning.

ROI

Bleaching experiment



Acquisition mode:
Type of imaging, 2D, 3D, time-lapse,
spectral scanning, simultaneous or
sequential imaging

Acquisition menus

Objectives

Auto mode. If you are unsure what
you are doing, start setting up your
imaging here.




Image compass
Lasers, detection windows, etc.

Live mode. Same as
"Capture Image" but on
continuous mode (speed
and resolution settings are
utilized here).

Fast live mode to check the
sample and find locations
with quick scanning (speed
and resolution settings are
NOT utilized here).

Detection configuration.
Three HyD detectors can be
turned ON-OFF, and their
parameters such as detection
range, gain and operating
mode can be modified.

Start to acquire 3D;
time-lapse;
or spectral scanning,
etc.


Capture an image



Auto mode (Dye Assistant)

Auto mode is a quick and easy way to create the settings for fluorophores and channels. Unexperienced users can start with the auto mode. Please note that auto mode might not work perfectly in special cases or has to be fine-tuned manually. Users shouldn't rely on auto modes, the safest way is to set up the imaging as manually as possible.

- Select your fluorophores
- Select between simultaneous or sequential imaging, line or frame






Manual mode, semi-manual mode (Image Compass)

Setting up the channels with the image compass is a love and hate thing. It was made user friendly and super-overcomplicated at the same time where you need to alternate between clicking buttons and dragging-dropping items. Watch out for high numbers of buttons, check-in marks, padlocks, menus and small signs. They all have a feature.

Sometimes you need to delete settings by dragging and dropping it to an empty space, sometimes to click on it and then click on the delete button.

Semi-manual mode:



1. Drag and drop your fluorophore to a detector
2. The excitation laser and detection window for the detector is selected automatically
3. Continue doing this with all your fluorophores starting from the blue range (HyD S 1 detector) till the red-infrared range (HyD S 3 detector)
4. To image sequentially, add a new "Setting" with the "+" sign and drag and drop a fluorophore to a detector in the new setting. Select "Stack, frame or line" mode.
5. Change the laser intensity by clicking on the laser line and changing the "intensity" (knob can be used)
6. Change the detector gain by clicking on the detector and changing the "gain" (knob can be used)


7. To delete a fluorophore from the settings, put your cursor on its emission curve and drag and drop it to an empty space. OBS ! If you want to only delete the laser line, click on it and push the trash button (top right corner)
8. Navigation between Setting 1 and 2. You can:
 - Select a setting by clicking on it
 - Uncheck (deactivate) a setting with the check-in box on the top left corner
 - Remove/delete a setting by the "X" button on the top right corner.

Manual mode:

Start with the semi-manual mode (by dragging and dropping a fluorophore to a detector) and modify the automatically generated channel (detection range, tuning the WLL onto a different laser wavelength, etc.).

OR

"Add Laser" (drag and drop this button to a desired laser wavelength) and **turn ON** a desired detector and **set up the detection window**, laser intensity and **detector gain**.



Get HELP and tutorial. Use it often.

It's encouraged to try out what all the different buttons and features do. Hover the mouse over a button to see the info about it.

If lost, don't forget to use the built-in HELP function of the software by clicking on the closest "i" button. This is one of the best features of LAS X.



Adjusting the correct intensity with the Detection configuration menu or the knobs

The best way to adjust the correct intensity of a channel is to use the Over-/Underexposure tool on the left side of the image window (first button). If Over-/Underexposure tool is activated we see the overexposed pixels in blue.

Over-/Underexposure tool

This channel is selected (clicked), we change the detector gain for this channel

Setting 1 | HyD S 1 | DAPI (H2O solved)

Setting 1 | HyD S 3 | ALEXA 568

Setting 2 | HyD S 2 | ALEXA 488

Preview: x = 512, y = 512 (768 kB)
Size: 369.82 μm x 369.82 μm

Ch1
Ch2
Ch3

3D

Overexposure (blue color); the histogram is clipped on the right side. To avoid clipping and to remove the blue, lower the "gain" or "laser power"

These channels are not selected (not clicked), we change the detector gain for the other channel

The overexposed pixels are out of the detector range, their intensity information is lost. We must avoid seeing overexposed (blue) areas. By changing (decreasing) the detector gain (voltage) or the laser intensity, we can move the overexposed areas into the range of the detector. For multichannel simultaneous imaging, be sure to click on the image of the channel you want to change, the other channel's detector is not changing meanwhile. Which one should we change, the laser or the



detector? There are no rules here, both have advantages and disadvantages. If we use high laser power with low gain we see a good quality image with low noise, but we can bleach the sample, so as always, COMPROMISE between quality and time/bleaching! Don't use a default laser or gain settings, always change them freely to get the best image without ruining the sample.

Note that if you change the pinhole or detection range, the signal is collected in a different Z/spectral range (intensity is changed), therefore new intensity adjustment is needed.

Acquisition menu

Optimize resolution (pixel size)

Set resolution (pixel size) manually. If you are unsure, use the optimize button. For thick 3D imaging use small, like 512x512 to save time/avoid bleaching.

Scanning speed (Hz): go up with the speed to decrease, go down to increase the quality. No golden rule, it is sample dependent, if unsure try the sample with different speeds, and use the best for the real acquisition. Again COMPROMISE between time and quality, decide what is more important.

Bidirecional scanning: Before starting a scan always test if X phase is correct. If unsure, don't use it, it can ruin you imaging.

Averaging: noise (random pixels) are removed with averaging process. Select the number of image to be used for averaging (again time vs. quality!) and choose between line and frame averaging. Use the check-in box to apply the settings to all sequences.

Pinhole: changes the thickness of the optical section. Click 1 AU button to get confocal imaging. If the signal is weak, open the pinhole more to collect photons from a thicker optical plane.

Format : 512 x 512

Speed : 400

Bidirectional X : OFF

Zoom Factor : 1.00

Zoom In : OFF

Image Size : 581.25 μm \times 581.25 μm

Pixel Size : 1.14 μm \times 1.14 μm

Optical Section : 3.121 μm

Pixel Dwell Time : 3.16 μs **Frame Rate :** 0.387/s

Line Average : 1

Line Accu : 1

Frame Average : 1

Frame Accu : 1

Rotation : 0.00

Pinhole


Unit : AU Airy 1

Pinhole : 1.00

Emission λ [nm] : 580



Z-stacking (3D imaging)



1. Go to Live or Fast Live mode, use the focus knob to find the end of the sample and push Begin to select the starting position of the sample where the first 2D stack should be made.


2. Focus to the other end of the sample and hit End to select the position where the last stack should be made.

3. Select the number of 2D stacks. Two ways of doing it, one can select Z-step Size and type in the interval between two stacks or select the Number of steps and type in the desired slice number. If unsure, go for the System Optimized button to not lose any information between two stacks.

4. Advanced optimization setup

5. Setup for linear laser or detector gain compensation when signal is getting lost in deep samples.

Deconvolution (Lightning)



With the Lightning module, on-the-fly deconvolution can be used to increase resolution and signal to noise ratio.


Use the slider to decide a balance between speed and resolution.

Strategy: Adaptive is always better, use this (numbers are based on and calculated for parts of the image with different signal to noise ratio).

Select the used mounting medium or select "custom" and type in manually the correct refractive index.



The image window



Right click on the image to Snapshot a view (for MIP or orthogonal view for example)



Saving, exporting

In LAS X projects can be created and saved. In a project, many images, datasets, snapshots can be saved and later re-used. Saved data will end up in .lif files that can be opened in LAS X, ImageJ, and Imaris.

Always save your image in RAW format (.lif) as it contains all the settings and information. If image file is needed, right click on a dataset and select Export.

The screenshot shows the 'Open projects' panel in LAS X. At the top, there are three buttons: 'New project', 'Open project', and 'Save all'. Below these is a toolbar with icons for 'New project', 'Open project', 'Save all', 'Save project', 'Export', 'Properties', and 'Show/Hide gallery'. The main area displays a project gallery with the following items:


Item	Size	Dimensions
Project	26.7 MB	
Image001	3.1 MB	xy
Image002	3.1 MB	xy
Image003	3.1 MB	xy
Image004	786 KB	xy
Image005	3.1 MB	xy
Image006	786 KB	xy
Image007	3.1 MB	xy
Image008	786 KB	xy
Series009	8.7 MB	xyz

Annotations with red arrows point from text boxes to specific UI elements:

- 'New project' button points to the 'New project' icon in the toolbar.
- 'Open project' button points to the 'Open project' icon in the toolbar.
- 'Save all' button points to the 'Save all' icon in the toolbar.
- 'Show/Hide gallery' button points to the 'Show/Hide gallery' icon in the toolbar.
- 'Save project' button points to the 'Save project' icon in the toolbar.
- 'Right click on an Image and select Properties to see all the imaging settings for the dataset. To re-use the settings for a new scanning, hit Apply settings in the pop-up window.' points to the 'Properties' icon in the toolbar.
- 'Right click on an Image and select Export if image files or video files are needed of the dataset' points to the 'Export' icon in the toolbar.


Shutting down the microscope

1. **Clean** after yourself, **put in** the smallest objective and **cover** the microscope stand with its dust cover.
2. **Use the logbook** and type in your imaging session's details.
3. In LAS X software **save** everything, and **turn OFF** all lasers



To turn OFF the lasers, use ONLY the ON-OFF sliders in the Laser Overview window (in the Acquire tab) or in the Laser Config window (in the Configuration tab).

4. **Close** down the LAS X software, **wait** till fully OFF
5. **Shut down** the computer, **wait** till fully OFF
6. **Turn OFF** the scanner box switches below the table:
 - i. Turn the laser interlock key to horizontal position (emission light OFF)
 - ii. Laser button OFF
 - iii. Power button OFF



6. **Wait 10-15 sec** so the instrument turns itself off. Don't turn off anything else that is not mentioned here.